# Ultrameter III<sup>™</sup> 9PTK Titration Kit

1111 1183

CONDUCTIVITY RESISTIVITY TDS ALKALINITY HARDNESS SATURATION INDEX ORP/FREE CHLORINE pH TEMPERATURE

> bluDock™ wireless data transfer available



... Since 1957

## Waterproof ULTRAMETER III

### Advanced Design • Superior Performance



- <u>FAST</u> one-touch measurements for conductivity, resistivity, TDS, ORP, free chlorine, pH and temperature
- <u>SIMPLE</u> user-intuitive prompts step you through alkalinity, hardness and LSI titrations
- <u>CONVENIENT</u> LSI Calculator accepts both grains and ppm hardness units
- EASY keypad calibration
- <u>STABLE</u> 4-electrode conductivity cell eliminates polarization, increasing accuracy with minimal maintenance
- POWERFUL microprocessor-based surface mount circuitry
- <u>UNIQUE</u> pH/ORP glass sensor is userreplaceable
- <u>VERSATILE</u> solution modes provide accuracy in diverse water quality applications
- <u>EFFICIENT</u> design features auto-off that minimizes power consumption — one 9V battery lasts up to 100 hours/5000 tests
- <u>ADVANCED</u> proprietary temperature compensation and TDS conversion algorithms
- <u>RELIABLE</u> non-volatile microprocessor provides data back-up, even when the battery is changed; this assures all calibrations and memory data are retained

## FAST READINGS • SIMPLE TO USE

### ACCURACY YOU CAN TRUST

Since 1957, the Myron L Company has designed and manufactured highly accurate and reliable analytical instruments for a wide variety of demanding applications. Thousands of professionals around the world rely on the performance of our instruments every day. Uses range from boiler water testing to ultrapure water control to medical instruments for artificial kidney machines. We are proud of the trust our handheld instruments and monitor/controllers have earned. Some of our instruments

have been in continuous use for over 40 years. Our latest digital handheld, the Ultrameter III 9P, is based on the tried and tested design of the Ultrameter II. The 9P adds the ability to perform in-cell alkalinity and hardness titrations either as stand-alone measurements or as part of an LSI calculation. Alkalinity, hardness, pH and temperature values can be manipulated in the LSI calculator to predict the effect of changes on water balance. Measurements are extremely accurate. And the 9P is simple to operate and maintain.

## 9 Parameters!

- Conductivity · Resistivity
- TDS Alkalinity Hardness
- Langelier Saturation Index
  with Calculator
- pH ORP & Free Chlorine
- Temperature

#### Accuracy You Can Trust

Like the Ultrameter II, the Ultrameter III delivers exceptional accuracy of ±1% of reading (±.1% at calibration point). This high level of accuracy is achieved through advanced four-electrode conductivity cell technology, a powerful microprocessor, proprietary circuit design and conversion algorithms for three of the most common solution types, as well as a unique pH/ORP sensor. With displayed values of up to 9999, the full four-digit LCD ensures resolution levels never before possible in such affordable instruments.

For accuracy in diverse applications, Ultrameter III allows you to select from 3 solution modes that model the behavior of the most common solution types: KCI, NaCI and 442<sup>™</sup> Natural Water. Additionally, the Ultrameter III features a user adjustable conductivity/ TDS conversion ratio for greater accuracy when measuring solutions not contained in the microprocessor.

Factory calibrated with NIST traceable solutions, each Ultrameter III may be supplied with both certification of traceability and NIST traceable solutions for definitive calibration.

#### Free Chlorine Readings

The Ultrameter III features the ability to convert Oxidation Reduction Potential measurements from mV to parts per million (ppm) free chlorine. The mV to ppm free chlorine conversion algorithm is based on a conversion curve developed through extensive testing performed by the Myron L Company. The accuracy of low-range ppm free chlorine readings has been increased by extrapolating from the published data. With this new feature Myron L instruments can measure a dynamic range of sanitizer concentrations that is wider than the range of a colorimetric test kit.

#### **Innovative Engineering**

The Ultrameter III is a prime example of how high-tech engineering can greatly simplify and streamline a task. Whether in the lab, industrial plant, or in a remote field location, electrometric measurements are accomplished in 3 easy steps:

- 1. Fill the cell cup
- 2. Push a parameter key
- 3. Take the reading

Temperature compensation and range selection are both rapid and automatic. User adjustable temperature compensation (up to 9.99%/°C) which can also be disabled for applications requiring non-compensated readings.

Digital results display instantly. Measurement is continuous, displaying real-time readings rather than a static snapshot of the solution.

#### Easy to Calibrate

calibrations All are quickly accomplished by pressing the  $\blacktriangle$  or  $\blacktriangledown$ keys to agree with our NIST traceable Standard Solution. When calibration is necessary, display prompts simplify pH calibration and make sure the correct buffer is being used. Plus, all parameters (excluding factoryset temperature) have an internal electronic setting that can be used for field calibration and as a check on pH/ ORP sensor life.

#### In-Cell Titration Functions

The Ultrameter III adds the ability to perform in-cell conductometric titrations that provides a convenient way to determine alkalinity, hardness and LSI in the field. This eliminates the need to collect and transport samples to another location for analysis. User intuitive display prompts guide you through titration procedures from start to finish. All required reagents and equipment are included in the 9P titration kit.

#### Water Balance Analysis

The Ultrameter III features both an LSI Calculator and an LSI Titration measurement mode. The Calculator allows you to perform what-if scenarios to predict how changes in solution parameters would affect the water balance of a system. The titration measurement function allows you to accurately calculate a saturation index value of a specific solution to determine whether the solution is balanced, scaling or corrosive.

#### Hardness Unit Conversion

The Hardness titration. LSI Titration and LSI Calculator functions allow you to set the hardness unit preference to either grains of hardness or ppm CaCO<sub>3</sub>.

#### **Wireless Benefits**

The optional bluDock<sup>™</sup> accessory package is an integrated data solution for your record keeping requirements, eliminating the need for additional hardware, wires and hassle, Because the user never touches the data, there is little opportunity for data tampering and human error. bluDock software has an easy to use interface with user intuitive functions for storing, sorting and exporting data.

#### More from Myron L

Myron L offers a full line of water quality instrumentation for every application. Call +1-760-438-2021 or visit our website for details: www.myronl.com

#### BENEFITS THAT SAVE YOU TIME & MONEY











Hardness/ Langelier Saturation Index Calculator allows you to determine water balance adjustments on the spot.

Easily transfer stored readings to Macintosh and PC platforms with the optional bluDock™ accessory package.

Ample memory provides increased flexibility to record and store 100 separate readings.

**Real Time Clock** with Date & Time Stamp allows you to maintain the integrity of each individual reading.

**In-Cell titrations** are simple to do. **Reagent solutions** are color-coded for easy recognition during procedure. Cell extender, pipette and replacement tips are included with kit.

The pH/ORP sensor chamber provides protection to a unique porous liquid-junction.

The large capacity KCI reservoir guarantees extended life.

A custom LCD helps simplify calibration and operation by using annunciators and prompts to indicate various conditions.

#### IP67/NEMA 6 rated

Ultrameter III is waterproof and buoyant and can be fully immersed to 3 feet/1 meter.



#### Features

- · Measures 9 Parameters: Conductivity, Resistivity, TDS, Alkalinity, Hardness, LSI, pH, ORP/Free Chlorine, Temperature
- LSI Calculator for hypothetical water balance calculations
- Wireless data transfer capability with bluDock<sup>™</sup> option
- · Autoranging delivers increased resolution across diverse applications
- · Adjustable Temperature Compensation and Cond/TDS conversion ratios for user-defined solutions
- Nonvolatile memory of up to 100 readings for stored data protection
- Date & time stamp makes record-keeping a snap
- · pH calibration prompts alert you when maintenance is required
- · Auto-off minimizes energy consumption
- · Low battery indicator

#### **Specifications**

\* 
 Sabic Innovative Plastics IP BV

Display	4 Digit Liquid Crystal Display				
Dimensions: LxWxH	196 x 68 x 64 mm/7.7 x 2.7 x 2.5 in.				
Weight	352 g/12.4 oz.				
Case/conductivity cell material	VALOX*				
Cell capacities (without cell extender)	pH/ORP: 1,2 ml/0.04 oz. Cond/TDS/Res: 5 ml/0.2 oz.				
Power	9V alkaline battery				
Battery life	>100 hours (5000 readings)				
Operating/storage temperature	0 – 55°C/32 – 132°F				
Protection ratings	IP67/NEMA6 Waterproof to 1 m/3 ft.				

#### Daramotore

	Conductivity	TDS	Resistivity	рН	ORP	Free Chlorine	Alkalinity Titration	Hardness Titration	LSI Titration	Temperature
Ranges	0–9999 µS/cm 10–200 mS/cm in 5 autoranges	0–9999 ppm 10–200 ppt in 5 autoranges	10 ΚΩ–30 ΜΩ	0–14 pH	±999 mV	0.00–10.00 ppm (5 <ph<9)< td=""><td>10-800 ppm</td><td>0-1710 ppm 0-100 grains</td><td>-10 to +10</td><td>0–71°C 32–160°F</td></ph<9)<>	10-800 ppm	0-1710 ppm 0-100 grains	-10 to +10	0–71°C 32–160°F
Resolution	0.01(<100 μS) 0.1(<1000 μS) 1(<10 mS) 0.01(<100 mS) 0.1(<200 mS)	0.01(<100 ppm) 0.1(<1000 ppm) 1(<10 ppt) 0.01(<100 ppt) 0.1(<200 ppt)	0.01(<100 KΩ) 0.1(<1000 KΩ) 0.1(>1 MΩ)	0.01 pH	1 mV	0.01 ppm	0.01(<100 ppm) 0.1(<1000 ppm) 1(<800 ppm)	0.01(<100 ppm) 0.1(<1000 ppm) 1(<1710 ppm)	0.01	0.1°C/F
Accuracy	±1% of reading	±1% of reading	±1% of reading	±0.01 pH*	±1 mV	± 2.5% of reading ppm	_	_	_	±0.1°C
Auto Temperature Compensation	0–71°C 32–160°F	0–71°C 32–160°F	0–71°C 32–160°F	0–71°C 32–160°F	-	_	_	_	-	-
Adjustable Temperature Compensation to 25°C	0–9.99%/°C	0–9.99%/°C	0–9.99%/°C	_	_	-	-	_	_	_
Conductivity/TDS Ratios Preprogrammed	KCI, 442**, NaCl	KCI, 442**, NaCl	_	_	_	_	_	_	_	_
Adjustable Conductivity/ TDS Ratio Factor	0.20-7.99	0.20-7.99	-	_	-	_	_	_	_	-

Certificates confirming the NIST traceability of an

Ultrameter III are available (must be specified when

Standard Solutions are necessary to maintain

accuracy and for periodic calibration of conductivity/

TDS parameters. All Standard Solutions are

NIST traceable for your complete confidence.

RECOMMENDED VALUES: KCI-7000 (conductivity),

Calibration Solutions are required for alkalinity,

hardness and LSI titrations. MODELS: Reagents

EDTA-HC available in 2 oz/59 ml and 1 qt/1 L.

Alkalinity Standard 100 and Hardness Standard

pH Buffers are necessary to maintain accuracy

and for periodic calibration of pH and ORP

parameters. Calibration with pH 7 Buffer is

especially important. All pH 4, 7, and 10 Buffers

H2SO4-1, H2SO4-3, NAOH, EDTA-LC and

200 available in 2 oz/59 ml, 1 gt/1 L, and 1

Limited Warranty

Hardness and Alkalinity Reagents and

placing instrument order). MODEL: MC9

442-3000 (TDS).

gal/3,8 L.

442 Natural Water Standard ™ Myron L Company

#### AHL Titration Kit

with hard protective case includes 9P; cell extender; titration plunger,  $100\mu$ L pipette; 12 tips; 3 buffers (pH 4, 7, and 10); pH/ORP Sensor Storage Solution; ORP Conditioner Solution; titration reagents H2SO4-1, HCI-3, NAOH, EDTA-LC, EDTA-HC;



Model 9PTK-BD Shown

Alkalinity Standard 100; Hardness Standard 200; standard solutions KCI-7000 and 442-3000. All bottles are 2 oz/59 ml. MODEL: 9PTK

9PTK AHL Titration Kit can be ordered with the bluDock Accessory Package, as well. KIT WITH BLUDOCK: 9PTK-BD

#### Accessories, Replacement Parts & Solutions

Soft protective case constructed of padded Nylon features a belt clip for hands-free mobility. MODELS: UCC (Blue); UCCDT (Desert Tan)

#### **Built on Trust**

Founded in 1957, Myron L Company is one of the world's leading manufacturers of water quality instruments. Because of our policy of continuous product improvement, changes in design and the specifications in this brochure are possible. You have our assurance any changes will be guided by our product philosophy: Accuracy, Reliability, Simplicity.



PRINTED IN U. S. A.

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE

2450 Impala Drive

are NIST traceable and are available in 2 oz/59 ml, 1 gt/1 L, and 1 gal/3,8 L.

pH Sensor Storage Solution available in 2 oz/59 ml, 1 gt/1 L, and 1 gal/3,8 L. MODELS: SS20Z, SSQ and SSG

Certificate of NIST traceability for pH Buffer or Conductivity Standard Solutions are available (must be specified when placing solution order). MODEL SC

#### **Replacement pH/ORP sensor**

user-replaceable; features a unique porous liquid-junction. MODEL: RPR

Myron L M'Pet 100µL pipette. MODEL: FVMP-100

Replacement Tips for Myron L M'PET micropipette. Bag of 48. MODEL: PTIP100-48

bluDock<sup>™</sup> Accessory Package includes bluDock<sup>™</sup>, Macintosh/PC application software for downloading data, dongle and printed instructions. MODEL: BLUDOCK

www.myronl.com

© MYRON L COMPANY 2011 UMITIDS 02-11

Tel: +1-760-438-2021 Fax: +1-800-869-7668 / +1-760-931-9189

Carlsbad, California 92010-7226 USA

